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## **INSTRUCTION FOR USE - READ CAREFULLY**

**Fresenius Kabi AG**  
61346 Bad Homburg  
Germany

**UK /AUS**

---

---

# **Glycine 1.5%**

Solution for Irrigation

---

---

<b>Properties</b>	Glycine is a naturally occurring amino acid.
<b>Composition</b>	1 liter contains: 2-amino acetic acid (Glycine) 15.0 g Water for injections in bulk, Ph. Eur.  pH value 4.5 – 6.5 Titration acidity 1.0 – 1.5 mmol NaOH/l Theoretical osmolarity 200 mosm/l
<b>Indications</b>	Intra- and postoperative irrigation during urological and gynecological procedures
<b>Contraindications</b>	Hyperhydration
<b>Precautions</b>	Caution is required in patients with pre-existing hyponatremia. Serum sodium levels must be corrected before use in this case.  Glycine 1.5% should be used with caution in patients with cardiac, renal or liver insufficiency.  Adequate precaution measures should be taken to avoid massive fluid absorption, e.g. balancing of the fluid in- and outflow (stopping the procedure if the difference exceeds 1 - 2 liters), limiting the operation time to 1 hour, and not exceeding the recommended maximum pressure for the intervention
<b>Side effects</b>	Fluid absorption may occur in operative transurethral and hysteroscopic procedures. In case of

---

---

---

---

intravascular absorption of large amounts of electrolyte-free irrigation solutions, hyperhydration and changes in the electrolyte balance (hyponatremia) may occur and consequently lead to acid base imbalance (metabolic acidosis). This condition is known as „TUR syndrome“. The associated symptoms include nausea, vomiting, increase or decrease in blood pressure, bradycardia, headache, restlessness, confusion, lethargy, chest pain, muscular twitches, dyspnoea, and poor urine output. In severe cases seizures, shock, pulmonary edema, cerebral oedema, cardiovascular collapse, and coma may occur. Fatal cases have been reported.

If the development of TUR syndrome is suspected the procedure should be terminated and adequate therapeutic measures initiated according to the individual condition of the patient.

Side effects associated with the properties and metabolism of glycine include visual disturbances (including transient blindness) and hyperammonaemia resulting in neurological disorders (including hyperammonaemic encephalopathy) in case of massive absorption.

**Interactions with other medicines**

None known.

**Dosage**

Dosage and quantity to be applied will result from the type of intervention.

**Description and connectivity**

The Careflex® Duobag is a polyolefin bag with a spiking port and the Care-Lock® port.

The carelock port is only compatible with the carelock transfer set.

The spiking port is compatible with all Fresenius Kabi transfer sets. Other spiking systems may be used and have to be tested for suitability by the user. Based on tests conducted with other transfer systems on the market the following sets are known to be not compatible with the spiking port of the Careflex® Duo

---

---

bag: Pump set for ARTHREX-System by Meise Medizintechnik (Ref: AR 6410), Pump Tube for fms solo-duo by Meise Medizintechnik (Ref: AR 7410), Tubing for Cystoscope or Resectoscope by Polfa-Lublin, Normoflo Irrigation Warming Set by Smiths Medical (Ref: IRI-600), AHTO Disposable Suction Irrigator Tube Set with Tip by Stryker (Ref: 250-070-620), Flocontrol Tube Set Pump SIL, Reusable by Stryker (Ref: Z1428-39), TUR Irrigation by Techno Med (Ref: TY30002SR), Fluidio Irrigation URO Set by The Surgical Company International B.V. (Ref: 670800-B), Fluidio Irrigation URO Spike Set by The Surgical Company International B.V. (Ref 680803-B), Single Use Tubing Set for Pump PV-5201 by Vimex Sp.Z.o.o. (Ref: PV-5201SUTS).

Use of non-compatible transfer systems may lead to fluid leaks or separation of the spike from the bag, during use, preventing the irrigation to continue.

An instruction how to spike the bag can be found at the end of this instruction for use.

**Presentation**

Polyolefin bags (Careflex® Duo)	5 x 2,000 ml
	4 x 3,000 ml
	2 x 5,000 ml

**Additional information**

Not for intravenous use or injection!

Do not use with bipolar surgical instruments.

The sterility of Glycine 1.5% solution for irrigation is guaranteed only until the container is opened. The solution must be used immediately after opening.

Use only if the solution is clear and free from particles and if the package undamaged.

For single use only.

Any unused portion of the solution must be discarded.

Re-use of medical single use equipment involves a risk of contamination of the equipment. This may

---

---

cause infection in the patient or the user, which may lead to illness or death.

Keep out of the sight and reach of children!  
Do not use Glycine 1.5% after the expiry date.

**Date of revision** June 2021

	CE marking of conformity with Medical Devices Directive 93/42/EEC (as amended) with the intervention of Notified Body 0123
	Indicates the presence of a sterile fluid path within the medical device in cases when other parts of the medical device, including the exterior, might not be supplied sterile.
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Indicates the medical device manufacturer.
	Indicates the need for the user to consult the instructions for use.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Indicates the date after which the medical device is not to be used.
	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Indicates a medical device that should not be used if the package has been damaged or opened.

---

---

	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	This product is not for intravenous use or injection.
	Indicates a medical device that needs protection from light sources.
	Indicates a medical device that needs to be protected from moisture.
	Indicates the temperature limits to which the medical device can be safely exposed. Store at room temperature.

---

---

---

---


## Spiking the bag

For spiking, the bag can be hung on an IV pole or alternatively placed on a flat surface. Administration of the content in hanging position.

Check product name, bag integrity and the expiry date. Inspect the container for damage or leakage and for particles. **Do not use if the container is damaged or if the solution contains particles.**




Identify the blue spiking port.

Break off the blue tamper-evident cover from the port.  
Do not touch the opened surface!

---

---

Close clamp of transfer set if applicable. Insert the spike **centered** and **straight** until the blue plastic collar of the port meets the shoulder of the spike.

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