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## Summary of the Study

### Aim of the study:

The aim of this study is to evaluate a potential deleterious effect of the FRESSENIUS KABI VOLUMAT AGILIA peristaltic pump for intravenous perfusion of human Red Blood Cell concentrates (RBC).

### Material and method:

The deleterious measured effect is the RBC's haemolysis when infused with VOLUMAT AGILIA pump and its associated tubing sets.

The protocol comprised two main steps:

1. First step : measurement of haemolysis and k<sup>+</sup> after one single pumping with the following variable parameters :
  - Age of the RBC's : 21 and 42 days
  - Flow rate : 200 ml/h and 1200 ml/h
  - Kind of pump : pump n°19800235 with standard membrane, and pump PSST 2-0110 with semi-rigid membrane (also called caterpillar).
2. Second step : multiple and successive pumping (10 times) is performed using a single outdated RBC (42 to 43 days) at 1200 ml/h. This trial was conducted on the two pumps.

The measured parameters are the haemolysis and the Extra-cellular potassium concentration.

### Results :

With a rate of haemolysis after one single pumping always below 0.1% and a very low release of potassium, RBCs do not undergo any significant deleterious effect by the pumping of the VOLUMAT AGILIA. The age of the RBCs (within a legal limit of 42 days), and the flow rate have not significant effect.

These results are confirmed by trials performed 10 times successively.

### Conclusion :

**The Volumat AGILIA Pump from Fresenius Kabi, used under normal supplier defined conditions do not alter the quality of the Human Red Blood Cell concentrates for transfusion.**

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