



CompoSeal Mobilea II

Portable tube sealer

Mobilea II meets the high demands of professional sealing and combines it with the benefits of modern battery technology.

High quality seals

Reliable operations

Quick and easy handling

Process control ensures leakage detection

CompoSeal Mobilea II

Reliable technology
for portable
and safe sealing

Intelligent power management

The CompoSeal Mobilea II utilises a modern Lithium Polymer battery for high sealing performance.

- Advanced Lithium Polymer battery provides a capacity of up to 1000 seals per full charge on standard blood bag tubes
- Designed for portable use: Low weight power pack of just 760g, robust metal housing, rubber pads for shock absorption
- High quality seals independent of the battery capacity
- Automatic adjustment of sealing time depending on tube thickness
- Approved for bedside sealing while donor is connected





CompoSeal Mobilea II is designed for safe and reliable tube sealing with the flexibility needed for mobile, in-house donation teams or processing departments.



CompoSeal Mobilea II Portable tube sealer

The portable light weight sealer with high battery capacity

Process efficiency to help you achieve more

Ordering information

For more information such as literature, technical details, working instructions and equipment, please contact your local sales representative.

CE marked, Class I according Directive 93/42/EEC, as amended by 2007/47/EC.

Electrical safety, IEC 60601-1

Cleared for sales in the United States



**FRESENIUS
KABI**

caring for life

Fresenius Kabi AG
61346 Bad Homburg
Germany

Phone: +49 61 72 608-0
Fax: +49 6172 608-5786

composeal@fresenius-kabi.com
www.fresenius-kabi.com