



**FRESENIUS
KABI**

caring for life



CompoCouch SALSA

The comfortable donation and apheresis chair

The new CompoCouch SALSA chairs feature modern ergonomics design with easy handling, mobility and stability for donors comfort and safety.

Compact and versatile for staff and donors

Hand control unit with initial positions button and manual locking

Stepless motorized adjustment from sitting to shock position

CompoCouch SALSA

CompoCouch SALSA offers a variety of models to meet your donors' needs for relaxed blood donations or apheresis procedures,

SALSA A1

Compact and versatile for staff and donor:

- Stepless motorized adjustment from sitting to shock position
- Foam armrests, folding upwards and rotatable
- Hand control unit with initial positions button, manual locking
- One motor for adjustment of complete seating



SALSA A2

Basic chair for donation and apheresis:

- Stepless motorized adjustment from sitting to shock position
- Foam armrests, folding upwards and rotatable
- Hand control unit with initial positions button, manual locking
- Two motors for adjustment of back and leg part



SALSA A3

Advanced chair for donation and apheresis:

- Stepless motorized adjustment from sitting, comfort and bed position to shock or Trendelenburg position
- Foam armrests, folding upwards, rotatable and tiltable
- Hand control unit with initial positions button, magnetic locking
- Three motors for adjustment of back, leg and seat part



SALSA A4

Premium chair for donation and apheresis:

- Stepless motorized adjustment from sitting, comfort and bed position to shock or Trendelenburg position
- Foam armrests, folding upwards, rotatable and tiltable
- Hand control unit with initial positions button, magnetic locking
- Four motors for adjustment of back, leg and seat part, plus height adjustment



CompoCouch SALSA, providing maximum support and comfort

CompoCouch SALSA chairs are easily adjustable and manageable for your donors' needs, with up to five positions available for SALSA chair models A3 and A4.



Sitting position



Comfort position



Bed position



Shock position



Trendelenburg shock position



HANDLING

All SALSA chairs have a hand control unit for comfortable adjustment



SAFETY

All SALSA chairs have safe working load of 200 kg



QUALITY

All SALSA chairs have a two-layer foam upholstery



COMFORT

All SALSA chairs have standard wedge shape arm rests and a neck roll



HYGIENIC

The framework of the SALSA chairs has homogenous surfaces to provide an easy and hygienic cleaning and with the paper roll holder, the neck and seat can always be freshly covered

CompoCouch SALSA Donation and apheresis chair

The choice is yours.



The CompoCouch SALSA is available in seven different colors.



Variety in accessories

- Transparent covers for foot part, foot step and armrest
- Castors with 10 cm Ø, individually lockable
- Foot rest
- Side table with tray

Manufacturer: LIKAMED GmbH, Raußmühlstrasse 7, 75031 Eppingen, Germany

CE marked. Class I according Directive 93/42/EEC.