



**FRESENIUS  
KABI**

caring for life

# CompoCouch SALSA

The comfortable donation and apheresis chair



## Technical Data

All SALSA models equipped with neck roll, wedge-shaped arm rests and paper roll holder. Neck roll and wedge-shaped armrest in upholstery color. Other upholstery colors for the SALSA models on demand.

Technical data	SALSA A1	SALSA A2	SALSA A3	SALSA A4
Size, length x width	199.5 x 80.8 cm	230.9 x 80.8 cm	229.3 x 80.8 cm	229.1 x 80.8 cm
Armrests, length x width x high rear/front	52.5 x 15.5 x 11.5/7.0 cm			
Mattress, length x width x height	210.0 x 53.5 x 7.0 cm			
Access height	64.0 cm	65.0 cm	63.5 cm	56.5 - 68.0 cm
Safe working load	200 kg			
Angle of incline, seat part	4°	9°	0°	0°
Power rating, primary	100 - 240 V, 50/60 Hz			
Electric safety	EN/IEC 60601-1			
Electromagnetic compatibility (EMC)	EN/IEC 60601-1-2			
Medical Device Class (93/42 EWG)	Class I			
Protection class under EN/IEC 60601-1	II			
Application part	Type B			
Protection type	IP X4			
CE classification 93/42 EWG (MDD)	CE			
Emission sound pressure level	< 65 dB (A)			
Weight	approx. 50 kg	approx. 55 kg	approx. 60 kg	approx. 60 kg
Environmental conditions:				
Temperature	10 - 40°C			
Humidity	30 - 75 % r.F.			
Air pressure	800 - 1060 hPa			

## Features

Features	SALSA A1	SALSA A2	SALSA A3	SALSA A4
Safe working load of 200 kg	■	■	■	■
Sitting, comfort and shock position	■	-	-	-
Sitting, laying and shock position	-	■	-	-
Sitting, comfort, bed, shock and Trendelenburg position	-	-	■	■
Stepless adjustment of positions	■	■	■	■
Hand switch for stepless adjustment	■	■	■	■
Folding upwards and rotatable wedge-shaped armrests	■	■	■	■
Tiltable wedge-shaped armrests	□	□	■	■
Paper roll holder	■	■	■	■
Individually lockable castors Ø 7.5 cm	■	■	■	■
Individually lockable castors Ø 10.0 cm	□	□	□	□
Comfortable soft upholstery	■	■	■	■
Neck roll in upholstery colour	■	■	■	■
Colour metall frame and cover Light Grey RAL 7047	■	■	■	■
Manually adjustable footrest	□	□	□	□
Colour footrest Dark Grey	■	■	■	■
Wedge-shaped armrest in upholstery color	■	■	■	■
transparent covers for leg part, footrest, armrest	□	□	□	□

■ Standard □ Optional

## Ordering Information

For more information such as literature, technical details and working procedures, please contact your local sales representative.

### Order number

9020941  
9020951  
9020901  
9020961  
9020911  
9020971  
9020921  
9020981

### SALSA models

SALSA A1, light grey  
SALSA A1, pastel blue  
SALSA A2, light grey  
SALSA A2, pastel blue  
SALSA A3, light grey  
SALSA A3, pastel blue  
SALSA A4, light grey  
SALSA A4, pastel blue

### Order number

LIKA-57-10-36-OPT  
LIKA-61-10-42-OPT  
LIKA-61-10-43-OPT  
LIKA-57-10-49-OPT  
LIKA-57-10-48-OPT  
LIKA-57-10-19-SET  
9028131

### Accessories

Protective cover for foot part SALSA, transparent  
Protective cover for foot rest SALSA, transparent  
Protective cover for arm rest SALSA, transparent  
Protective cover for low, wedge-shaped arm rest Salsa, transparent  
Protective cover for wedge-shaped arm rest Salsa, transparent  
Foot rest SALSA, with clamping, for backfitting  
Side table with tray

Manufacturer: LIKAMED GmbH Raußmühlstrasse 7, 75031 Eppingen Germany.

CE marked. Class I according Directive 93/42/EEC.



**FRESENIUS  
KABI**

caring for life

Fresenius Kabi AG  
61346 Bad Homburg  
Germany

Phone: +49 61 72 608-0  
Fax: +49 6172 608-5786

compocouch@fresenius-kabi.com  
www.fresenius-kabi.com