

Data sheet **Base A**

DATA CONCENTRATOR



The Base A is a Data Concentrator and supplies mains power for up to 8 Module DPS and/or Module DPS Visio and/or Module MVP.

The Base A enables the connection of Orchestra Infusion Workstation to a Patient Data Management System via a unique RS 232 cable.



**FRESENIUS
KABI**

caring for life

BASE A

According to the number of intravenous infusions needed, at any combination, up to 8 Modules (Orchestra Module DPS or Module DPS Visio or Module MVP) can be fitted onto the Orchestra Base A.

A single mains cord and a single communication cable help reduce installation constraints.

TECHNICAL SPECIFICATIONS	
MAINS SUPPLY	From 95 to 240 V - 50/60Hz.
BATTERY	Characteristics: 6V 1.1/1.3 Ah - sealed lead acid battery rechargeable. Battery life: 2 hours minimum.
ALARMS, PREALARMS AND SAFETY	The Base A gives a synthetic view of the operation of connected Modules.
COMPLIANCE, STANDARDS	Complies with the 93/42 EEC European Medical Device Directive: CE 0459. Complies with EN 60601-1, EN 60601-1-2 and EN 60601-2-24 requirements.
WATERPROOFNESS	Protection against ingress of liquid IP34 (CEI 529).
ELECTRICAL SAFETY	Protection against leakage current: type CF equipment. Protection against electric shocks: class II equipment.
COMMUNICATION INTERFACE	Compatible with major PDMS software.
RS232 INTERFACE	One unique cable for all connected Modules.
DIMENSIONS / WEIGHT	H 31,5 x W 5,5 x D 20 cm / 2 kg.

Orchestra

Infusion Workstation composed of Bases and Infusion Modules adapted to every clinical practice.

Module DPS
Modular Syringe Pump



Module DPS Visio
Modular Syringe Pump



Module MVP
Modular Volumetric Pump



Base A
Data concentrator



Base Intensive
Unique 2 channel relay function, Customised drug library, Infusion history management.



Base Primea
Simultaneous control of both hypnotic and analgesic effects in TCI.



Orchestra Mobile stand
Whatever may be the architecture of your department, Orchestra Mobile Stand provides the choice of customised accessories.



Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for current situation.

Manufacturer:

Fresenius Vial S. A. S.
Le Grand Chemin
38590 Brézins
France
Tel.: +33 (0)4 76 67 10 10
Fax: +33 (0)4 76 67 11 34

A Company of Fresenius Kabi

Disposables

A wide range for use with Module DPS and Module MVP. Catalogue is available from our Customer Service.